


**DNA ladder HMW****Cat #BB-D30 (500µl)****Cat #BB-D31 (500µl×2)**

This ready to load high molecular weight (HMW) DNA ladder with eleven bands (12kb, 9kb, 7kb, 6kb, 5kb, 4kb, 3kb, 2kb, 1.5kb, 1kb, 0.5kb) is premixed with Xylene Cyanol loading dye.

**Storage: -20°C****Note:** 5µl loading is recommended.**HMW**  
0.8% TAE Agarose Gel