



BioBharati LifeScience Pvt. Ltd.

Developing Best Products for Quality Research of National Significance

EN-35, First Floor, Salt Lake, Sector-V, Kolkata, INDIA, PIN: 700091

Registered Office: 79, Purbachal Kalitala Road, Kolkata, INDIA, PIN: 700078

E-mail: biobharati@gmail.com, Web: www.biobharati.com

Product Name

Cat. No.

Pack Size

DNA Ladder LMW^{Plus}

BB-D40

50 µg (500 µL)

BB-D41


2 X 50 µg (1 mL)

Concentration: 0.1µg/µL

Storage: -20°C

This ready to load low molecular weight DNA ladder contains thirteen bands (3000bp, 2000bp, 1500bp, 1000bp, 900bp, 800bp, 700bp, 600bp, 500bp, 400bp, 300bp, 200bp, and 100bp). This DNA ladder premixed with loading dye (Bromophenol blue and Xylene Cyanol) can be used for direct loading on gel (~5 µL/ lane).

Note: 5µL (i.e., 500ng) loading is recommended.



LMW ladder visualized by EtBr staining on a 1% TAE agarose gel. Mass value are for 500 ng/lane.